

## **CEM Thermometer Operating Instructions**

- For the first use or when inserting new batteries wait between 10 minutes for the warm-up of the thermometer and when inserting the new batteries
- When the thermometer is initially turned on the LCD screen will read “LO”
- Once it is pointed to someone’s forehead from a distance of between 2 inches and 5.9 inches it will read the temperature. Make sure there is no hair, perspiration, cosmetic or cap covering the forehead prior to taking the temperature.
- Make sure the thermometer is in “Body” mode
- When the batteries are low there will be a battery symbol that will display on the LCD screen
- When the screen reads “LO” it may be due to a temperature of less than 89.6 degrees Fahrenheit in “Body” mode or less than 32 degrees Fahrenheit in “Surface” mode